

SIGH BULLETIN


JULY 2006
For Theatre and Sterile Services Managers

SIGH acquire Osborn and Simmons

Pictured: Examples from the Osborn and Simmons range of high-quality, ophthalmic instruments now available from SIGH



Foldable Lens Folding Forceps



Cross Action Lens Insertion Forceps





Direct Insertion Forceps




Corneal Section Scissors - Left



Corneal Section Scissors - Right




Vannas Scissors



Acquisition will provide customers with outstanding levels of service

The recent acquisition of Osborn and Simmons has made SIGH one of the UK's leading specialist providers of high quality, **ophthalmic surgical instruments**. The Osborn and Simmons range includes ophthalmic forceps, scissors, needleholders, choppers, manipulators, pushers and spatulas. Catalogues and prices are available and the whole range is featured on SIGH's website at www.sighltd.co.uk



OSBORN & SIMMONS

- Disposable Cannulae
- Phaco Instruments
- Capsulorhexis Forceps
- Micro Forceps
- Lens Dialers
- Nucleus Manipulators
- Nucleus Choppers
- Nucleus Rotators
- Lens Introducing
- Lens Folding & Insertion Forceps
- Capsulotomy Scissors
- Corneal Scissors
- Suture Tying Forceps
- Markers
- Tying Forceps
- Needleholders
- Iris Retractors
- Sterilizing Trays

Colour coded locking tags



Shown here actual size

A practical and effective marking system of brightly coloured tags that quickly identify or mark instruments, carts, trays, cabinets and other equipment.

Features: Available in 12 distinctive colours • Non-sterile Autoclavable • Tamper resistant • Latex Free Durable and economical • Quick and easy to attach Versatile - adapt to your own specific coding needs • Brightly coloured tags replace hard to see safety pins and sutures


Click!

Just £18.00 per 100 tags



Re-usable Bipolar & Diathermy forceps

We supply re-usable Bipolar and Diathermy forceps for all specialities. We also supply leads for both types of forceps. Contact us at the number below for a full colour brochure



DIA-021 Diathermy Forceps

BP-200-ISF Bipolar Forceps

Surgical instrument cleaning brushes

General instrument cleaning brush
7.25" with stainless steel bristles for use on rasps, files, reamers, burrs, needle holders, jaws, vascular serrations and stubborn stains.

General cleaning brush
9" utility brush for use on sinks, autoclaves, ultrasonic cleaners, bowls, basins, I.V. poles, commodes, wheelchairs etc.



General instrument cleaning brushes
single and double-ended, 7" and 7.25" with nylon bristles for use on scissors, hemostats, needle holders, dressing and tissue forceps.

Upright instrument cleaning brush
5.25" with fine stainless steel bristles for use on rasps, files, reamers, burrs, needle holders, jaws, vascular serrations and stubborn stains.

Upright instrument cleaning brush
5.25" with nylon bristle for use on scissors, hemostats, needle holders, dressing and tissue forceps.

Repair service for surgical instruments

- Re-insulation of Bipolar, Diathermy and Laparoscopic Instruments, Thudicums, Retractors and Skids.
- Sharpening of scissors, knives, osteotomes and chisels.
- Realignment of forceps, serrating and insertion of new T/C inserts in needleholders.
- Repair of Sphygs, Flowmeters and suction equipment.
- Bipolar and Diathermy **Healthcheck** for hospitals who must check their forceps after 50 autoclavings. SIGH collect and check the insulation and integrity of the forceps and cables. Instruments which pass are returned free of charge. Any instruments needing re-insulation are taken care of and failed cables are replaced and all returned on a weekly turnaround basis. We use BOWA cables and replace forceps with our own ISO registered product.

SIGH

Surgical Instrument Group Holdings Limited 89a Gloucester Road Croydon CR0 2DN
Telephone 020 8683 1103 Lo Call 0845 310 1103 Fax 020 8683 1105 Email info@sighltd.com



www.sighltd.com