


Colour coded locking tags



Shown here actual size

A practical and effective marking system of brightly coloured tags that quickly identify or mark instruments, carts, trays, cabinets and other equipment.

Features: Available in 12 distinctive colours • Non-sterile
 Autoclavable • Tamper resistant • Latex Free
 Durable and economical • Quick and easy to attach
 Versatile - adapt to your own specific coding needs • Brightly coloured tags replace hard to see safety pins and sutures

Click!

Just £18.00 per 100 tags



High quality fibre optic light cables

Universal fibre optic cables with adaptors for all light sources and instruments. Up to 35 possible combinations.

- 6', 7.5' and 9' lengths
- 1,2,3,4,5 and 7mm diameters

- Special lengths available
- Colour range: Grey, Pink, Purple, Blue, Black, Red, Green and Gold
- All light source and instrument ends
- Autoclavable


■ Prices for a complete cable, light source and instrument adaptor assembly start at just £225 plus carriage and VAT

Bipolar and Monopolar Leads

Manufactured by **BOWA Electronic GmbH**, who use the technology they have developed to make the electrical cabling for BMW and Volkswagen cars. These cables, in their distinctive grey colouring, are sealed units and so avoid fluid ingress which is the chief cause of cable failure and deterioration.



Diathermy Quivers



Detailed leaflet available on request

- 95-010-10 Julian Taylor 20cm round base with drainage holes
- 95-010-20 20cm X 6cm flat base
- 95-010-30 Laparoscopic Quiver 38cm X 6cm

For more information on the products featured in this Bulletin or for help with ordering, please call

0845 310 1103

All calls charged at local rates

Surgical instrument cleaning brushes

Our comprehensive range of cleaning equipment is designed for use with all kinds of surgical instruments including cannulae, tubes, laparoscopic instruments, orthopaedic reamers and the Laryngeal Mask.

Examples from the wide range of cleaning brushes available

General instrument cleaning brush 7.25" with stainless steel bristles for use on rasps, files reamers, burrs, needle holders, jaws, vascular serrations and stubborn stains.

General cleaning brush 9" utility brush for use on sinks, autoclaves, ultrasonic cleaners, bowls, basins, I.V. poles, commodes, wheelchairs etc.




General instrument cleaning brushes single and double-ended, 7" and 7.25" with nylon bristles for use on scissors, hemostats, needle holders, dressing and tissue forceps.

Upright instrument cleaning brush 5.25" with fine stainless steel bristles for use on rasps, files reamers, burrs, needle holders, jaws, vascular serrations and stubborn stains.

Upright instrument cleaning brush 5.25" with nylon bristle for use on scissors, hemostats, needle holders, dressing and tissue forceps.

Re-usable Bipolar & Diathermy forceps

We supply re-usable Bipolar and Diathermy forceps for all specialities. We also supply leads for both types of forceps. Contact us at the number below for a full colour brochure



DIA-021 Diathermy Forceps

BP-200-ISF Bipolar Forceps

Specialist surgical instruments from around the world

SIGH, together with its daughter company, Caterham Surgical Supplies, are U.K. distributors for some of the leading ranges of specialist surgical instruments from around the World.

Contact SIGH for a wide range of general surgical instruments for most specialities and for **MicroFrance** laparoscopic instruments.

SIGH also provide scissors, dressing, artery and dissecting forceps, rongeurs and ENT micro instruments.

Contact Caterham Surgical Supplies for:

Lawton general instruments

Lepine gynae instruments

Bowa leads and finger switches

You can reach SIGH and Caterham Surgical Supplies on the following LOCAL telephone number: **0845 310 1103**. All calls charged at local rates.



SIGH

Surgical Instrument Group Holdings Limited 89a Gloucester Road Croydon CR0 2DN

Telephone 020 8683 1103 Lo Call 0845 310 1103 Fax 020 8683 1105 Email info@sighltd.com



www.sighltd.com